

# PerKit™ HRP–Oligo Conjugation Kit (CM53407) User Reference Guide

## Contents

Important Notes & Contact Information .....	2
Kit Configuration and Components .....	3
Safety Information .....	4
Labeling Chemistry.....	4
Support .....	5
Protocol.....	6
1. Lab Instrumentation Needed.....	6
2. Preparation of Oligo Sample for Labeling.....	6
3. Convert Amine-Oligo to Disulfide-Oligo (Step 1).....	7
4. Convert Disulfide-Oligo to Thiol-Oligo (Step 2).....	8
5. HRP Conjugation with Thiol-Oligo (Step 3).....	9
6. Purification to Remove Excess HRP .....	9
Other Considerations.....	11
1. Concentration Determination for Oligo (Unlabeled).....	11
2. Concentration Determination for HRP–Oligo Conjugate.....	11
3. MW Calculation.....	11
4. Recommended Storage Conditions .....	12
5. Characterization of HRP–Oligo by HPLC .....	12
6. Sample Submission for HPLC Analysis .....	12
Appendix: Typical Kit Performance Data (HPLC analysis, CellMosaic).....	13



CellMosaic, Inc.  
10A Roessler Road  
Woburn, MA 01801  
USA

Phone: 781-463-0002  
Fax: 781-998-4694  
Email: [info@cellmosaic.com](mailto:info@cellmosaic.com)  
Website: [www.cellmosaic.com](http://www.cellmosaic.com)

---

## Important Notes & Contact Information

### READ BEFORE USING ANY RESOURCES PROVIDED HEREIN

The information provided in this document and the methods included in this package are for information purposes only. CellMosaic provides no warranty of performance or suitability for the purpose described herein. The performance of labeling using this kit may be affected by many different variables, including but not limited to purity and complexity of the HRP, differences in preparation techniques, operator abilities, and environmental conditions.

Sample data are provided for illustration and example purposes only and represent a small dataset used to verify kit performance in the CellMosaic laboratory. Information about the chemicals and reagents used in the kit are provided as necessary.

### For Research Use Only. Not for Use in Diagnostic Procedures.

The information in this document is subject to change without notice. CellMosaic assumes no responsibility for any errors that may appear in this document. In no event shall CellMosaic be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple, or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

### NOTICE TO PURCHASER: LIMITED LICENSE

The purchase of this product includes a limited, non-transferable license to use this product to practice the labeling methods using the reagents solely for the purchaser's research activities. The license granted herein is personal to the original purchaser and may not be transferred to any other party outside the purchaser's company. No other right or license is conveyed or granted either expressly, by implication or by estoppel, to resell or repackage this product. Further information can be obtained by contacting:

Director of Licensing  
c/o CellMosaic, Inc.  
10-A Roessler Road, Woburn, MA 01801.  
Phone: 781-463-0002  
Fax: 781-998-4694  
E-mail: [info@cellmosaic.com](mailto:info@cellmosaic.com)

## Kit Configuration and Components


CellMosaic® has designed this PerKit™ HRP–Oligo conjugation kit to prepare single-labeled HRP–oligo conjugates with linkers and conjugation chemistry that retains the activity of HRP. The same configuration and chemistry have been tested for HRP-oligos with different oligo sequences at a few diagnostic companies since 2011.

This kit is an upgraded and improved replacement kit with slightly different components and workflow as the first-generation kit (Cat# CM53401). This kit has multiple configurations for various lengths of oligos and can be used to conjugate 5 nano-mole (nmol) of oligo. **Table 1** provides the catalog numbers for various kit configurations and **Table 2** lists the kit components.

**Table 1:** Configurations of the PerKit™ HRP–Oligo conjugation kit (CM53407)

Configuration	No. of Reactions	Catalog No.
1. Labeling for short oligo (5–30 bases)	1	CM53407.1x1
	3	CM53407.1x3
2. Labeling for medium oligo (31–59 bases)	1	CM53407.2x1
	3	CM53407.2x3
3. Labeling for long oligo (≥60 bases)	1	CM53407.3x1
	3	CM53407.3x3

**How to use this protocol:** The protocol in this user manual is written for three different configurations: 5–30 bases oligo (CM53407.1), 31–59 bases oligo (CM53407.2), and ≥60 bases oligo (CM53407.3). Steps are common to all configurations. However, some of the kit components are specific to each configuration.

 Please follow the specific instructions for individual configurations.

### Requirement for amine-oligo (≥5 bases):

1. Total amount: 5 nano-mole oligo content as measured by UV.
2. HPLC purified and lyophilized oligo with amine-oligo content ≥90% by C18 HPLC.
3. ≥5- bases (oligos longer than 60 bases may result in lower loading and high impurities).

**Note:** You can purchase HPLC purified C6 amine-modified oligo from standard oligo manufacturers, up to 90 bases for this kit preparation. Oligos tested at CellMosaic for this kit were from Integrated DNA technologies (<https://www.idtdna.com/>) or TriLink Biotechnologies (<https://www.trilinkbiotech.com/>). We have not tested oligos made by other oligo manufacturers. It is highly recommended that customers analyze and quantify oligos prior to using them.

**Disclaimer:** CellMosaic has no financial interest in IDT or TriLink.

**Table 2:** Components and storage temperatures for PerKit™ HRP–Oligo conjugation kit. All kits share the same components except for the number of tubes provided for Reagent A and the type of Filter Devices.



- Upon receipt, please remove **Box 1** and store in a freezer at or below -20°C.
- Store **Box 2** in a refrigerator at 2–8°C.

Box No. (Storage T)	Name	Cat#	Kit Cat#	Quantity (x1)	Quantity (x3)
<b>Box 1</b> (≤-20°C)	C6 Maleimide Activated HRP (red label)	CM53214	Any	1 unit	3 units
	Reagent A (white label)	CM12002.1	CM53407.1	1 unit	3 units
			CM53407.2	2 units	6 units
Reagent B (blue label)	CM13001	Any	1 unit	3 units	
<b>Box 2</b> (2–8°C)	Solution A (green label)	CM01003	Any	0.3 mL	1 mL
	Buffer A (orange label)	CM02001		4 mL	12 mL
	Buffer B (indigo label)	CM02005		4 mL	12 mL
	Buffer C (equilibration and washing buffer, lime label)	CM02003		4 mL	10 mL
	Buffer D (elution buffer: 50 mM Tris buffer, pH 8.0, 1M NaCl, navy label)	CM02004		1.5 mL	5 mL
	Filter Devices (Only one part# will be supplied, 3 or 9 units)	CM03CD003A	CM53407.1	3	9
		CM03CD010A	CM53407.2 CM53407.3		
	Collection Tubes for Filter	CM03CT0	Any	6	18
	Column Q	CM03SC5		1	3
Collection Tubes for Column Q	CM03CT6	2		6	
User Supplied Material	Amine Oligo (≥5 bases)	N/A	NOT PROVIDED (1–5 nano-mole) for each reaction)		

## Safety Information


Warning: Some of the chemicals used can be potentially hazardous and can cause injury or illness.

Please read and understand the Safety Data Sheets (SDS) available at [www.CellMosaic.com](http://www.CellMosaic.com) before you store, handle, or use any of the materials.

## Labeling Chemistry

The kit is designed to work with amine-modified oligos. The user supplies the amine-modified oligos, which are readily available from many commercial oligo suppliers. Using the kit components, the user converts the amine-oligo to a thiol-oligo, followed by reaction of the thiol-oligo with activated HRP to generate the HRP-oligo conjugates. The total length of the linkage between the HRP and the amine-oligo

is 29 atoms excluding the linker used for amine modification of the oligo. The Q-column purification step typically provides the resulting HRP-oligo at greater than 80% purity.



**Scheme 1.** Synthetic route to HRP–oligo conjugate.

An HPLC purified and lyophilized oligo is required for this conjugation. If the content of the amine oligo is >90%, the purity of the final conjugate will generally be >90% pure with an average of 1 oligo per HRP (majority are single labeled). This purity is sufficient for ELISA, and a trace amount of impurities will negligibly affect the signal intensity. For oligos longer than 60 bases or oligos modified with non-phosphate backbone or bases, the purity of the amine oligo may be lower. In addition, if there is a solubility issue in aqueous buffer, the labeling reaction may be inefficient, resulting in low loading and high impurities. For complete removal of oligo impurities, gel filtration chromatographic purification is recommended.

Alternatively, you can use PerKit™ HRP oligo conjugation with a commercial-ready disulfide-labeled oligo (Cat#: CM53402x1 or CM53402x3 ) to improve the loading.

Key features of this HRP-oligo conjugation kit:


- High quality maleimide-activated HRP for the conjugation: >99% purity and >200 units/mg protein activity.
- A single purification affords 85-90% single-labeled HRP-oligo conjugates if the quality of the oligo is at par.
- Preparation can be done in a day.
- Stable linkage with optimal spacer to ensure no interference between oligo binding and HRP activation.
- All reagents included, from preparation to purification.
- Options to choose tailored services at CellMosaic after conjugation:
  - You can choose to send your conjugates to CellMosaic for HPLC analysis of the sample or complete removal of trace oligo impurities and unreacted HRP.

## Support

A customer can request a recommendation for the conjugation if the oligo has some special features or solubility issues. CellMosaic also provides fee-based support services to customers who need help analyzing the final conjugates by HPLC and further purification to remove trace oligo impurities.

## Protocol

There are a total of three reactions in this protocol. Keep in mind that thiol-labeled oligo should be purified quickly once the reaction is done and used immediately. HRP-oligo conjugation is generally done within 4 h, but you can take longer or leave it at 2–8°C overnight after the reaction. The Q-column purification will take less than 10 minutes.



**Scheme 2.** Schematic diagram of the workflow for preparing HRP–oligo conjugates.

### 1. Lab Instrumentation Needed

- Vortex mixer, centrifuge (preferably refrigerated, 14,000 g capable), mini-centrifuge
- Pipettes and tips
- Timer
- Incubator or shaker set at 37°C and RT
- Support stand, lab frame, or any support rod for desalting column
- Flask
- Personal protection equipment (lab coat, safety glasses, and chemical-resistant nitrile gloves)
- UV spectrophotometer (optional)


### 2. Preparation of Oligo Sample for Labeling

Items needed: Amine-Oligo (user supplied), 1 Filter Device (CM03CD003A or CM03CD010A), 2 Collection Tubes (CM03CT0), Buffer A (CM02001, orange label), clean centrifuge tubes (not provided in the kit).

Total amount of oligo used for the conjugation is **2–5 nano-mole**.

**A1.** Insert a new Filter Device into one of the provided Collection Tubes. Perform the step based on the following conditions.

- ✓ If your oligo is supplied as a lyophilized solid, dissolve the oligo in deionized water to a concentration of 50  $\mu\text{M}$  and then transfer up to 100  $\mu\text{L}$  to the Filter Device. Add Buffer A (orange label) to make up the total volume to 500  $\mu\text{L}$  and cap it.
- ✓ If your oligo is supplied as liquid, transfer up to **5 nano-mole** to the Filter Device directly. Add Buffer A to make up the total volume to 500  $\mu\text{L}$  and cap it.



**A2.** Place the capped Filter Device into the centrifuge rotor, aligning the cap strap toward the center of the rotor; counterbalance with a similar device.

**A3.** Spin the Filter Device at 14,000  $\times g$  (preferably cooled to 4°C) for 8 to 15 minutes to concentrate to < 100  $\mu\text{L}$ . Spin time will depend on which Filter Devices are supplied in the kit.

Catalog No# (1 rxn, 3 rxn)	Oligo Length	Spin Time	Typical Leftover Vol.
CM53407.1(x1, x3)	5-30 bases	15	80 $\mu\text{L}$
CM53407.2(x1, x3) CM53407.3(x1, x3)	$\geq 60$ bases	8	35 $\mu\text{L}$
Typical leftover volume is obtained using an Eppendorf 5417R and centrifuge at 4°C.			

**A4.** Remove the assembled device from the centrifuge and separate the Filter Device from the collection tube. Transfer the filtrate from the collection tube to a clean centrifuge tube (not provided). **Save the filtrate until the experiments are done.**

**A5.** Insert the Filter Device back into the collection tube. Add 400  $\mu\text{L}$  of Buffer A to the Filter Device. Spin the device at 14,000  $\times g$  for 8 to 15 minutes to concentrate to < 100  $\mu\text{L}$ . Remove the assembled device from the centrifuge. Transfer the filtrate from the collection tube to a clean centrifuge tube (not provided). **Save the filtrate until the experiments are done.**

### 3. Convert Amine-Oligo to Disulfide-Oligo (Step 1)

Items needed: Amine-Oligo from Step A5, 1 or 2 tubes of Reagent A (CM12002.1, white), Buffer A (CM02001, orange label), 1 Filter Device (CM03CD003A or CM03CD010A), 2 Collection Tubes (CM03CT0).

**B1.** Briefly spin the tube containing Reagent A (white label). Make sure you can see a tiny droplet at the bottom of the tube.

**B2.** Transfer the Amine-Oligo sample from the Filter Device from **Step A5** to the Reagent A tube by pipetting (preferably using a sterilized 100  $\mu$ L pipette tip with filter). Make sure the pipette tip reaches the bottom of the Filter Device.

**B3.** Wash the Filter Device twice with 10  $\mu$ L of Buffer A and transfer the wash to the same Reagent A tube from **Step A5** (**Note: Wash = add buffer, aspirate with pipette 2-3 times**).

**B4.** Vortex the reaction mixture for 10 seconds to mix and then spin down. **If two tubes of Reagent A are supplied in Box 1 (Cat# CM53407.3 for labeling  $\geq$ 60 bases oligo), transfer the mixture from the 1<sup>st</sup> Reagent A tube to the 2<sup>nd</sup> Reagent A tube. Repeat the mixing by vortexing the mixture for 10 seconds and then spin down.**

**B5.** Incubate at 37°C for 1 h.



Start Time: \_\_\_\_\_ End Time: \_\_\_\_\_

**Tip for mixing:** You can use a nutator, shaker, vortex, or incubator shaker for mixing. If you are using end to end nutating, make sure the tube from **step B4** is securely capped. If you don't have any of this equipment, you can let the tube sit on the bench with manual mixing by pipetting every 10 minutes.

**B6.** Remove the Tube from **Step B5** from the incubator and spin to ensure no liquid is in the cap.

**B7.** Insert a new Filter Device into a new provided Collection Tube. Transfer the reaction mixture from **Step B6** into the Filter Device directly. Wash the centrifuge tube twice with  $\sim$ 200  $\mu$ L of Buffer A, transfer the solution to the Filter Device (total volume should be  $\leq$ 500  $\mu$ L), and cap it. Spin the Filter Device at 14,000 x g for 8 to 15 minutes to concentrate to < 100  $\mu$ L (refer to **Step A3** for spin time).

**B8.** Remove the assembled device from the centrifuge and separate the **Filter Device** from the collection tube. Transfer the filtrate from the collection tube to a clean centrifuge tube (not provided). **Save the filtrate until the experiments are done.**

**B9.** Insert the **Filter Device** back to the collection tube. Add 400-450  $\mu$ L of **Buffer A** to make up the total volume to 500  $\mu$ L. Spin the device at 14,000 x g for 8 to 15 minutes to concentrate to < 100  $\mu$ L (refer to **Step A3** for spin time).

**B10.** Repeat **steps B8** and **B9** one time.

#### 4. Convert Disulfide-Oligo to Thiol-Oligo (Step 2)

Items needed: Disulfide-Oligo from Step B10, Reagent B (CM13001, blue), Buffer A (CM02001, orange label), Buffer B (CM02005, indigo), 1 Filter Device (CM03CD003A or CM03CD010A), 2 Collection Tubes (CM03CT0).

**C1.** Centrifuge the tube containing Reagent B (blue label) to spin down the solid. Make sure you can see a small amount of solid at the bottom of the tube.

**C2.** Transfer the Disulfide-Oligo sample from the Filter Device from **Step B10** to the Reagent B tube by pipetting (preferably using a sterilized 100  $\mu$ L pipette tip with filter). Make sure the pipette tip reaches the bottom of the Filter Device.



**C3.** Wash the Filter Device twice with 10  $\mu$ L of Buffer A and transfer the wash to the Reagent B tube from **Step C2 (Note: Wash = add buffer, aspirate with pipette 2-3 times).**

**C4.** Incubate at 37°C for 30 min.



Start Time: \_\_\_\_\_ End Time: \_\_\_\_\_



The following steps are to be performed without any break. Reduced thiols tend to oxidize quickly. Work quickly through **Steps C5-C7 and D1-D3.**

Work quickly

**C5.** Insert a new Filter Device into a new provided Collection Tube. Transfer the reaction mixture from **Step C4** into the Filter Device directly. Wash the centrifuge tube twice with  $\sim$ 200  $\mu$ L Buffer B (indigo label), transfer the solution to the Filter Device (total volume should be  $\leq$ 500  $\mu$ L), and cap it. Spin the Filter Device at 14,000 x g for 8 to 15 minutes to concentrate to < 100  $\mu$ L (refer to **Step A3** for spin time).

**C6.** Remove the assembled device from the centrifuge and separate the **Filter Device** from the collection tube. Transfer the filtrate from the Collection Tube to a clean centrifuge tube (not provided). **Save the filtrate until the experiments are done.**

**C7.** Insert the **Filter Device** back into the Collection Tube. Add 400-450  $\mu$ L of **Buffer B** to make up the total volume to 500  $\mu$ L. Spin the device at 14,000 x g for 8 to 15 minutes to concentrate to < 100  $\mu$ L (refer to **Step A3** for spin time).

## 5. HRP Conjugation with Thiol-Oligo (Step 3)

Items needed: Activated HRP (CM53211, red), Activated Oligo from Step C7, Solution A (CM01003, green).

Work quickly

**D1.** Add 125  $\mu$ L of Solution A (green label) to the tube containing Activated HRP (red label). Vortex for 30 seconds to 1 minute to dissolve the HRP.

**D2.** Transfer the Thiol-Oligo sample from the Filter Device from **Step C7** to the Activated HRP tube by pipetting (preferably using a sterilized 100  $\mu$ L pipette tip with filter). Make sure the pipette tip reaches the bottom of the Filter Device.

**D3.** Wash the Filter Device twice with 20  $\mu$ L of Buffer B and transfer the wash to the Activated HRP tube from **Step D2 (Note: Wash = add buffer, aspirate with pipette 2-3 times).**

**D4.** Incubate at 25°C or room temperature **in the dark (wrap the tube in aluminum foil or place in a closed incubator)** for at least 4 h, preferably overnight.



Start Time: \_\_\_\_\_ End Time: \_\_\_\_\_

## 6. Purification to Remove Excess HRP

Items needed: HRP-Oligo Reaction Mixture from Step D4, Buffer C (CM02003, Lime), 1 Column Q (CM03SC5), Buffer D (CM02004, navy), 2 Collection Tubes (CM03CT6).

**E1. Sample Preparation:** When **Step D4** is complete, add 175  $\mu\text{L}$  of Buffer C (equilibration and washing buffer, lime label). Pipette the solution up and down in the tube three times to mix.

**E2. Column Equilibration:** Add 400  $\mu\text{L}$  of Buffer C to a Column Q and centrifuge at 2000 x g for 2 minutes. Discard the flow-through. Repeat this step one more time.

**E3. Sample Application:** Add the sample solution from **Step E1** to the equilibrated Column Q at a maximum volume of 400  $\mu\text{L}$ . Centrifuge at 2000 x g for 5 minutes. Set aside the flow-through.

**Tip:** Pay attention to the color of Column Q; if the conjugation reaction is successful, the top of Column Q will be light brown. **If the top of Column Q is colorless, please store your remaining buffers/solvents/Column Q inside a 4°C refrigerator and contact CellMosaic.**

**E4. Washing:** Add 400  $\mu\text{L}$  of Buffer C to Column Q and centrifuge at 2000 x g for 2 minutes. Repeat two more times. Make sure the final wash is colorless. Discard the flow-through.

**E5. Elution:** Place Column Q on a clean Collection Tube. Add 400  $\mu\text{L}$  of Buffer D (elution buffer, navy label) to the column. Centrifuge at 2000 x g for 2 minutes. Collect the flow-through and label it as elution (**Elution 1**).

**Tip for elution volume:** **Elution 1** contains the majority of your purified HRP-oligo conjugate. However, conjugating longer oligos may require a greater volume to elute. Pay attention to the color of Column Q; if it is still brown-colored, repeat **Step E5** to elute and collect more fractions.

### HRP-Oligo is Ready for Your Experiment

**Tip:** The approximate concentration of **Elution 1** is 6  $\mu\text{M}$  in 400  $\mu\text{L}$  of 50 mM Tris buffer, pH 8.0, 1 M NaCl. Elution can be diluted (for lower concentration) or a desalting column (not included in kit) can be used to exchange the buffer.

---

## Other Considerations

### 1. Concentration Determination for Oligo (Unlabeled)

The accuracy of the amount of oligo is important for achieving optimized oligo labeling in this protocol. The oligo manufacturer will usually supply the oligo in lyophilized form with the amount measured prior to lyophilization. This quantitation is usually sufficient. For accuracy, you can re-measure the concentration of the oligo after re-suspending the oligo in deionized water in **Step A1**.

$$\text{Concentration (M) of oligo before dilution} = \frac{(A_{260}) \times DF}{L \times \epsilon_o}$$

**A<sub>260</sub>**: UV absorbance of oligo at 260 nm.

**L**: UV cell path length (cm) - if you are using a 1 cm UV cell, dilute the oligo in water to a concentration of 1 to 2 μM.

**DF**: dilution factor.

**ε<sub>o</sub>**: extinction coefficient of oligo.

### 2. Concentration Determination for HRP–Oligo Conjugate

To determine the concentration of your conjugate, dilute your HRP–oligo from **Step E5** with 1x PBS buffer. Measure the UV absorbance of the HRP–oligo at 403 nm (A<sub>403</sub>) using a UV spectrometer and calculate the concentration using the following formula:

$$\text{Concentration (}\mu\text{M) of the diluted sample} = \frac{A_{403} \times 10}{L \times 1.02}$$

$$\text{Concentration (mg/mL) of the diluted sample} = \frac{A_{403} \times 10 \times (43000 + \text{Mw(oligo)})}{L \times 1.02}$$

**L**: UV cell path length (cm) - if you are using a 1 cm UV cell, you can dilute HRP–oligo 5 to 10-times to get a good reading

**Mw(oligo)**: molecular weight of oligo.

### 3. MW Calculation

Calculation of the MW of the conjugate:

$$\text{Mw(conjugate)} = \text{Mw(oligo)} + 43000$$

#### 4. Recommended Storage Conditions

For long-term storage, HRP-oligo conjugates can be lyophilized and stored as lyophilized powder at -80°C for 1 year.

#### 5. Characterization of HRP–Oligo by HPLC

HRP–oligo conjugate can easily be characterized by size-exclusion chromatography (SEC) and anion exchange chromatography (AEX) HPLC. SEC separates the molecules by apparent molecular weight (MW) or size in aqueous solution. The larger the MW of the conjugate, the earlier it elutes. By comparing the SEC profile of the unlabeled HRP, unlabeled oligo, and conjugate, you can check whether a conjugate formed, the heterogeneity of the labeling, and the percentage of unlabeled oligo and HRP. AEX separates the molecules based on the net surface negative charge (total charge). After conjugation, the total negative charge of the conjugates may be different and result in separation. It may be difficult to obtain good data with AEX HPLC without optimization. We recommend customers use SEC HPLC as the main method for conjugate analysis.

CellMosaic offers two SEC standards ([Product #: CM92004](#) and [CM92005](#)) for our customers to use with any SEC column. The CM92004 product sheet contains all of the information and methodology needed to run an SEC HPLC analysis.

If you do not have access to an HPLC facility, you can send your sample to CellMosaic for analysis.

#### 6. Sample Submission for HPLC Analysis

If you are submitting samples to CellMosaic for HPLC analysis, please follow these instructions:


- 1) Go online: <https://www.cellmosaic.com/hplc-analysis/>, select SEC HPLC Analysis ([Product# AS0023](#)) and/or AEX HPLC Analysis ([Product# AS0026](#)), choose the quantity (number of samples - bulk discounts are available for multiple samples), and submit the order. Alternatively, you can email [info@cellmosaic.com](mailto:info@cellmosaic.com) for a quote and to place the order.
- 2) Dilute your un-conjugated oligo to 10-20 µM in PBS buffer, and then transfer 50 µL of the diluted solution to a 0.5 mL micro-centrifuge tube. Label the vial properly.
- 3) Dilute your conjugate to 0.1–0.5 mg/mL in PBS buffer, and then transfer 50 µL of the diluted solution to a 0.5 mL micro-centrifuge tube. Label the vial properly.
- 4) Ship your samples with a cold pack for overnight delivery.

## Appendix: Typical Kit Performance Data (HPLC analysis, CellMosaic)


**Oligo information:** HPLC-purified 18 bases oligo with 5' NH<sub>2</sub>-C6 modification

**Kit Lot number:** S453.S8.0801C (Box 1)

**Figure 1:** Size-exclusion HPLC analysis of the C6 Maleimide Activated HRP (Detect at 260 nm).



**Figure 2:** Size-exclusion HPLC analysis of the oligo (Detect at 260 nm).



**Figure 3:** Size-exclusion HPLC analysis of purified conjugate (detect at 260 and 400 nm).

